
# Wasserstrahlpumpe – Water Jet Filter Pump

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## Operating Manual

Please read carefully before use,  
and follow all operating and safety  
instructions!

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## Contents


	Page
<b>1. Safety Instructions</b> .....	<b>4</b>
<b>2. Operating Limitations</b> .....	<b>6</b>
<b>3. Components</b> .....	<b>7</b>
<b>4. Application</b> .....	<b>7</b>
Possible applications of the filter pump .....	7
<b>5. Connection to water supply</b> .....	<b>8</b>
Threaded pipe connection.....	8
Hose connection .....	8
Performance check of the non-return valve .....	9
Connection to vacuum line.....	9
<b>6. Cleaning</b> .....	<b>10</b>
<b>7. Replacing non-return valve</b> .....	<b>11</b>
<b>8. Technical Data</b> .....	<b>12</b>
Table Ultimate pressure.....	14
Table Suction capacity.....	14
Table Water consumption .....	14
<b>9. Warranty</b> .....	<b>15</b>
<b>10. Ordering Data</b> .....	<b>16</b>

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## Safety Instructions

### Please read the following carefully!

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this pipette to consult and establish appropriate safety and health

practices and determine the applicability of regulatory limitations prior to use.

1. Every user must read and understand this operating manual prior to using the instrument and observe these instructions during use.
2. Follow general instructions for hazard prevention and safety instructions; e.g. wear protective clothing, goggles and gloves. Observe special instructions for handling evacuated vessels, and if necessary use a protective shield.
3. Only use evacuated vessels which are definitely made for work under vacuum (e.g. desiccators and filter flasks).

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\*DVGW = German Association of Gas and Engineering

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
## Operating Limitations

1. The user should always ascertain whether the filter pump is suited for the particular application.
2. Some chemicals or their fumes may damage this filter pump (swelling, embrittlement). Do not use with chemicals which attack NBR, FKM, PP and PTFE.
3. Avoid exposure to temperatures exceeding 80 °C. Precision moulded parts of plastic material tend to distort at these temperatures.
4. The supply water pressure must be between 3 and 6 bar. At a lower pressure, the specified performance cannot be attained. A higher pressure may damage the pump.
5. The function of the non-return valve may be impaired through deposits or impurities (⇒ page 9).

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## Application


### Possible applications of the filter pump:

- Acceleration of filtering processes (vacuum filtration)
- Gentle distillation or sublimation (pressure reduction lowers the boiling or sublimation point).
- Drawing off liquids.
- Evacuation, e.g. of desiccators.

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
## Connection to water supply:

### 1. Threaded pipe connection

- Screw 3/4" sleeve nut tightly on water tap.
- A thread adapter for 1/2" with gasket, NBR is included. For 3/8" and M 22 x 1, adapters are also available. (⇒ page 16)
- Screw the required thread adapter tightly on the water tap. Then connect the filter pump. Do not use force.

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

### 2. Hose connection

- Unscrew the water inlet nozzle (with sleeve nut) from the water pump.
- Screw the water inlet nozzle (with hose nipple) tightly on the pump.
- Slide the water supply hose on the tube nipple till it fits tightly. Secure the water supply hose with a clip.

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
### 3. Performance check of the non-return valve

- Connect the filter pump to the water supply pipe (without vacuum line), open water inlet – the pump sucks.
- Close the water outlet by pressing your finger tip on it.
- No water should run out of the hose nipple now (vacuum arm).
- If there is, however, water leaking, renew the non-return valve (⇒ page 11).

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### 4. Connection to vacuum line

- Unscrew the hose nipple (vacuum arm) before sliding on the vacuum tubing.
- Slide the vacuum tubing onto the loose hose nipple until it fits tightly.
- Re-fit the hose nipple on the filter pump.

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#### Important:

Do not slide on the vacuum tubing as long as the filter pump is connected to the water tap.

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

### Removal of internal sediments

- Unscrew the water inlet nozzle (thread or hose nipple).
- Unscrew the hose nipple for vacuum connection.

### Hinweis:

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### Note:

Do not dismantle the filter pump any further!

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
- Immerse all parts into diluted hydrochloric acid (1:1) or a common deliming agent.
- Allow time for reaction, then rinse the parts thoroughly with water.
- Re-assemble all parts in the reverse order. Screw tight by hand. Do not use force.

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## Replacing non-return valve

In case of water leakage during testing procedure of the non-return valve, the valve disc and PTFE sealing ring have to be replaced:

- Remove the hose nipple (vacuum arm) and O-ring.
- Pull out the PTFE sealing ring with a suitable tool (e.g. angled forceps).
- Remove the valve disc.
- Insert a new valve disc. **Make sure that it is correctly placed in front of the valve opening.**
- Slide in a new PTFE sealing ring right to the stop.

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For new non-return valve please refer to page 16, "Ordering data".

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water

### Connection to water supply:

R 3/4"; with thread adapter R 1/2";  
with hose nipple: 10 - 12 mm external diameter

### Connection to suction line:

hose nipple, 6 - 9 mm external diam.,  
with screw cap GL 14

### Length of the filter pump:

approx. 210 mm (with 3/4" connection)

### Weight:

approx. 33 g (with 3/4" connection)

### Material:

Flat packings: NBR (nitrile butadien rubber)  
O-rings: FKM (Fluoro elastomer) (An O-ring made of  
NBR fixed in the housing.)  
Valve reed: FKM (Fluoro elastomer)  
Sealing ring: PTFE  
All other parts: PP

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– continued –

### Water consumption

approx. 190 litres/h at 3.5 bar (absolute pressure)

### Suction capacity

400 l/h ( $\pm$  50 l/h) air against atmospheric pressure at a supply pressure of 3.5 bar and a water temperature of 12 °C.


### Ultimate pressure

at 12 °C (water temperature):  $\leq$  16 mbar


### Water supply pressure

to attain the ultimate pressure: 3 - 6 bar (absolute)

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
## **Warranty**

We shall not be liable for the consequences of improper handling, use, servicing, operation or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original parts have been used.

Subject to technical modification without notice. Errors excepted.

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### Items supplied:

#### Water jet filter pump, complete

with sleeve nut R 3/4",  
thread adapter from R 3/4" to R 1/2",  
water inlet nozzle (hose nipple 10 - 12 mm Ø)  
and vacuum connection (hose nipple 6 - 9 mm Ø),  
with valve reed of FKM (Fluoro elastomer)

<b>Cat. No.</b>	1596 00
-----------------	---------

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### Hose nipple

for vacuum connection, 6 - 9 mm Ø,  
with sleeve nut and sealing, FKM

<b>Cat. No.</b>	1596 69
-----------------	---------

### Non-return valve

with valve reed of FKM,  
complete

<b>Cat. No.</b>	1596 68
-----------------	---------