
**Dynalab Corp.**


175 Humboldt Street, Suite 300  
Rochester, NY 14610  
Tel: 585-334-2060 800-828-6595  
Fax: 585-334-0241  
[www.dynalab.com](http://www.dynalab.com)

# Stuart Melting Point Apparatus SMP10

Instructions for use  
Manuel d'utilisation  
Istruzioni per l'uso  
Instrucciones de funcionamiento  
Bedienungsanweisung

# English

**Figure 1**

**Thank you for purchasing this Stuart product. To get the best performance from the equipment, and for your own safety, please read these instructions carefully before use.**

This equipment is designed to operate under the following conditions:-

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuation not exceeding 10%
- ❖ Overvoltage category II IEC664
- ❖ Pollution degree 2 IEC664

If the equipment is not used in the manner described in this manual the protection provided by the equipment may be impaired.

# Electrical Installation

 **THIS EQUIPMENT MUST BE EARTHED**

**Before connection please ensure that the line supply corresponds to that stated on the rating label. Power consumption is 70W.**

There is an IEC socket at the rear of the right hand side of the instrument for connection to the mains supply, (see figure 1). The unit is supplied with two mains leads fitted with IEC plugs for connection to the instrument. One has a U.K. 3 pin plug and the other has a 2 pin "Shuko" plug for connection to the mains supply. Choose the lead appropriate for your electrical installation and discard the other.

Should neither lead be suitable you should obtain a moulded lead locally. If this is not possible, take the lead with the U.K. plug and replace the plug with a suitable alternative. See the enclosed instruction sheet for advice on how to carry out this procedure.

Should the mains lead require replacement a cable of 1mm<sup>2</sup> of harmonised code H05W-F connected to an IEC 320 plug should be used. N.B. The UK mains lead is protected by a 10A fuse mounted in the plug top.

The mains lead should be connected to the instrument BEFORE connection to the mains supply.

**IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN**

# General Description

The Stuart SMP10 has been designed for maximum safety and ease of use. The temperature is selected, measured and displayed digitally ensuring accuracy and avoiding the need for a glass thermometer. Two samples can be tested simultaneously. They are viewed via a magnifying lens with clear observation aided by built in illumination. Extendible feet allow the unit to be operated at the optimum viewing angle. Full access to the block allows easy cleaning and maintenance.

## Before Using your SMP10 for the first time

Once the unit is unpacked, select the appropriate choice of self adhesive instruction label from the multi-language set provided. Peel off the backing and stick the label onto the left side of the SMP10 fascia. Apply the label carefully using a flat object to smooth away air bubbles.

# Principle of Operation

The SMP10 melting point apparatus has been designed to give both quick and accurate results as well as maximum convenience in use.

Melting point samples are placed in a glass capillary tube which is placed in the aluminium block inside the sample chamber. This block is heated and the sample observed through the magnifying lens until the melt occurs. The melting point temperature may then be easily read from the large LED display.

In order to avoid the necessity to continually watch the sample the SMP10 is equipped with a "plateau" function. This allows a temperature to be set a few degrees below the expected melting point. The SMP10 will then heat to this temperature very rapidly (20°C per minute) and hold it until the operator is ready to begin measuring. It can safely be left to do this unattended.

Once measuring is started the SMP10 will heat slowly (2°C per minute) from the plateau temperature until the melt occurs. This slow rate of heating allows very accurate melting points to be obtained.

### **Before Use**

Place the unit on a firm, level surface and extend the feet at the rear of the instrument so that the unit is at a convenient viewing angle.

Connect to the mains electricity supply and switch on by means of the mains on/off switch. The actual block temperature will now be shown on the three digit LED display and the illumination in the sample chamber will come on.

### **Measuring a Melting Point**

N.B. If the unit has been used recently the block may be too hot for your sample. If this is the case press the stop button and allow the unit to cool before proceeding.

1. Prepare sample by placing a small amount in the end of a glass capillary tube.
2. Decide on a suitable plateau temperature. This should be approximately 10°C below the expected melting point of the sample.
3. Check that all 3 function lights are extinguished. If not, press the stop button.
4. Press and hold the set button (the plateau light will flash). The display will now show the current plateau temperature. When first switched on this will be 50°C but note that if the unit has recently been used the value may have been re-set by the previous user.
5. The desired plateau temperature may now be set using the arrow keys to scroll the display up or down as required.
6. Release the set button. The new plateau temperature is now set and all function lights will go out (The plateau setting can be checked at any time during operation by pressing and holding the set button). This action will not interfere with the operation of the unit).
7. Insert the tube into the side of the heating block via the holes provided. For convenience this can be done from either side of the block. Look down the magnifier and position the tube so that the sample can be observed clearly.  
**Caution: The heating block may be HOT.**
8. Press the start key. The unit will quickly heat up to the plateau temperature (the heating light will be illuminated - see figure 1)
9. Once the plateau temperature has been reached the plateau light will be illuminated as well as the heating light - see figure 1. N.B. At low plateau temperatures there may be some overshoot. Wait until the plateau light comes on before proceeding in order to ensure that the temperature has stabilised.
10. Press the start button again. The block will begin to heat at the pre-set ramp rate of 2°C per minute (The plateau light will go out and the ramping and heating lights will both be illuminated - see figure 1.)
11. Observe the sample until the melt occurs and record the temperature from the digital display.
12. After the melt has occurred press the stop button. All function lights will go out and the unit will cool to ambient temperature.
13. Pressing the start button again will cause the unit to return to the plateau temperature instead of ambient. The heating light will come on even if the temperature is above the plateau and the unit is, in fact, cooling.

# Cleaning & Servicing

**WARNING:** Ensure the unit is disconnected from the mains electricity supply and allowed to cool before attempting any cleaning or servicing.

## Cleaning

Periodically clean the instrument using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents.

Broken melting point tubes may be removed from the block simply by pushing them through with a piece of wire or similar tool.

For access to the block for more thorough cleaning the sample chamber cover must be removed.

1. Remove the 3 retaining screws and lift off the cover.
2. Loosen the screw retaining the metal plate but do not remove completely.
3. Push the plate to one side.
4. The glass window in the top of the block may now be removed.

## Servicing

Any repairs or replacement of parts MUST only be undertaken by suitably qualified personnel.

There are two replaceable fuses on the rear of the IEC mains connector. The cause of any fuse failure should be investigated and corrected.

Proceed as follows:

- 1) Ensure the SMP10 is disconnected from the mains supply.
- 2) Unscrew the four screws on the base of the unit and remove the base plate. Do not disconnect or pull on the base earth bonding wire.

- 3) Two fuses are mounted on the rear of the IEC connector (both 'Live' and 'Neutral' mains supply leads are fused). These fuses can be removed and replaced by hand. ALWAYS replace the fuse by ones of the correct values shown below:

For 230V, 50Hz - T1A

For 120V, 60Hz - T1.6A

- 4) Refit the base and four securing screws.
- 5) It is advisable to perform a 'PAT' (Portable Appliance Test) electrical safety test, or equivalent, before the SMP10 is put back into operation.

For a comprehensive list of parts and a Service Manual required by service engineers conducting internal repairs, please contact the Technical Service Department of Bibby Sterilin Ltd. quoting both the model and serial number.

Only spare parts supplied or specified by Bibby Sterilin Ltd. or its agents should be used. Fitting of non-approved parts may affect the performance and safety features designed into the instrument.

If in any doubt, please contact the Technical Dept. of Bibby Sterilin or the point of sale. Bibby Sterilin Ltd. Stone, Staffordshire ST15 0SA, United Kingdom.  
Tel: +44 (0) 1785 812121  
Fax: +44 (0) 1785 813748  
e-mail [equipment@bibby-sterilin.com](mailto:equipment@bibby-sterilin.com)

# Warranty

Bibby Sterilin Ltd. warrants this instrument to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of two (2) years. In the event of a justified claim Bibby Sterilin will replace any defective component or replace the unit free of charge.

This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by incorrect installation, adaptation, modification, fitting of non-approved parts or repair by unauthorised personnel.


**Bibby Sterilin Ltd.**

Stone, Staffordshire, ST15 0SA  
United Kingdom  
Tel: +44 (0) 1785 812121  
Fax: +44 (0) 1785 813748  
E-mail: [equipment@bibby-sterilin.com](mailto:equipment@bibby-sterilin.com)  
Internet: [www.bibby-sterilin.com](http://www.bibby-sterilin.com)

# Technical Specification

|                        |                                                 |
|------------------------|-------------------------------------------------|
| Number of samples      | Two                                             |
| Temperature range      | Ambient to 300°C                                |
| Temperature resolution | 1°C                                             |
| Temperature accuracy   | ±1°C of displayed temperature                   |
| Display                | Three digit LED                                 |
| Sensor                 | PT100                                           |
| Ramp rates             | 20°C / minute to plateau<br>2°C /minute to melt |
| Dimensions (WxDxH)     | 160 x 220 x 170 mm                              |
| Net weight             | 1.8 Kg                                          |

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Bibby Sterilin Ltd. Stone, Staffordshire ST15 0SA, United Kingdom.

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## Notes

These products meet the relevant EC harmonised standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in their vicinity will meet these standards

and we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, please contact the Technical Service Department of Bibby Sterilin Ltd.

## ***Declaration of Conformity***

---

### **Model: SMP10 Melting Point Apparatus**

This product complies with the requirements of the EU Directives listed below:

**89/336/EEC** Electromagnetic Compatibility Directive, amended by 92/31/EEC and 93/68/EEC.

**73/23/EEC** Low Voltage Directive amended by 93/68/EEC

Compliance with the requirements of these Directives is claimed by meeting the following standards:

**EN61326-1:1997** Electrical equipment for measurement control and laboratory use  
Emissions and Immunity

**EN61000-3-2:1995** Mains harmonics

**EN61000-3-3:1995** Voltage flicker

**EN61010-1:1993** Safety requirements for electrical equipment for measurement, control and laboratory use. Plus part2-010:1995 particular requirements for laboratory equipment for heating of materials

The technical documentation required to demonstrate that the products meet the requirements of these Standards is available for inspection on request.

Signed:  (Mr C. R. Baggaley)

Date: 14/2/00

Authority: Technical Manager

for

Bibby Sterilin Ltd - Stone - Staffs - ST15 OSA - UK  
Tel: +44 (0) 1785 812121 - Fax +44 (0) 1785 813748

# **INSPECTION REPORT**

**MODEL** SMP10

## **ELECTRICAL SAFETY**

1. Earth continuity
2. Insulation
3. Flash test


## **FUNCTIONAL**

1. Indicators
2. Temperature control
3. Visual acceptance

**QUALITY CONTROL INSPECTOR**

### **Dynalab Corp.**

350 Commerce Drive  
Rochester, NY 14623

Tel: 585-334-2060 800-828-6595

Fax: 585-334-0241

[www.dynalabcorp.com](http://www.dynalabcorp.com)