

Fits up to 72 pegs  
**High Impact PS Drying Rack**



Dynalon<sup>®</sup>  
Labware



# Interchangeable pegs Wall mount or benchtop



259184

Back view of  
259184 in  
159184-0001



259184  
drain outlet



- Injection-molded, high-impact, light blue polystyrene
- Pegs hold cylinders, bottles, beakers and containers etc. and snap in or out easily when desired
- Comes with 72 large 4 in. long (10.2cm) pegs, expansion hanging hooks, draining tubes and fastening stoppers for coupling with additional racks
- Interchangeable pegs for use with test tubes are available separately
- Bag of 11 Pegs made of high impact polystyrene
- Dimensions: 17W x 24 in. H (43.2 x 61cm) with 3 in. (7.6cm) ledge

## Description

259184	Rack Drying 72 Place, PS
259194	Rack Bag of 11 Pegs for 5918-0000
159184-0001	Drying Rack Stand, Acrylic

