




Polystyrene Container Vials with Snap-on Caps

Part #	Description	Material	Gram Wt	OFC	A	B	C	Bulk Pack	Gross Pack
7003001	Vial, PS 3Dr, Bulk	PS	3.8+/-0.5gr.	10cc+/-10%	1.900+/-0.015	.790+/-0.015	.670+/-0.015	2,250	
8003001	Cap, White, LDPE 3Dr, Bulk	LDPE	.6+/-0.1gr.	n/a	.255+/-0.015	.910+/-0.020	.765+/-0.015	14,000	
426364-03	Container w Snap Cap, PS 3Dr								144
7005011	Vial, PS 5Dr, Bulk	PS	4.0+/-1gr.	18cc+/-10%	1.973+/-0.015	.970+/-0.015	.850+/-0.015	1,400	
8003010	Cap, White, LDPE 5Dr, Bulk	LDPE	1.0+/-0.2gr.	n/a	.255+/-0.015	1.080+/-0.020	.930+/-0.015	10,000	
426364-05	Container w Snap Cap, PS 5Dr								144
7007010	Vial, PS 7Dr, Bulk	PS	6.8+/-1gr.	25cc+/-10%	1.970+/-0.015	1.150+/-0.020	1.020+/-0.015	1,000	
8003020	Cap, White, LDPE 7/10Dr, Bulk	LDPE	1.2+/-0.3gr.	n/a	.275+/-0.015	1.280+/-0.020	1.110+/-0.015	7,500	
426364-07	Container w Snap Cap, PS 7Dr								144
7010020	Vial, PS 10Dr, Bulk	PS	10.5+/-1gr.	38cc+/-10%	3.210+/-0.015	1.132+/-0.015	1.000+/-0.015	600	
426364-10	Container w Snap Cap, PS 10Dr								144
7015001	Vial, PS 15Dr, Bulk	PS	9.5+/-1gr.	52cc+/-10%	2.435+/-0.015	1.430+/-0.015	1.315+/-0.015	500	
8003050	Cap, White, LDPE 15Dr, Bulk	LDPE	1.6+/-0.3gr.	n/a	1.560+/-0.020	1.558+/-0.015	1.389+/-0.015	5,000	
426364-15	Container w Snap Cap, PS 15Dr								144
7016020	Vial, PS 20Dr, Bulk	PS	13.0+/-1gr.	69cc+/-10%	2.480+/-0.015	1.660+/-0.015	1.510+/-0.015	400	
8003060	Cap, White, LDPE 20Dr, Bulk	LDPE	2.7+/-0.2gr.	n/a	.300+/-0.015	1.800+/-0.015	1.630+/-0.015	3,000	
426364-20	Container w Snap Cap, PS 20Dr								144
7030001	Vial, PS 30Dr, Bulk	PS	16.5+/-0.2gr.	113cc+/-5%	2.923+/-0.015	1.925+/-0.015	1.790+/-0.015	230	
8003070	Cap, White, LDPE 30/40Dr, Bulk	LDPE	3.4+/-0.3gr.	n/a	.305+/-0.015	2.065+/-0.015	1.900+/-0.015	2,500	
426364-30	Container w Snap Cap, PS 30Dr								72
7040001	Vial, PS 40Dr, Bulk	PS	26.0+/-1gr.	145cc+/-5%	3.842+/-0.015	1.935+/-0.015	1.780+/-0.015	175	
426364-40	Container w Snap Cap, PS 40Dr								72