

Unlimited classroom uses
Student Project Tanks



Dynalon
Labware


PRICE
DROP



**Aquariums, terrariums or critter chambers,
this series is ready to meet your needs!**

All accessories are included.

Shatter-resistant polycarbonate Safety in the classroom!



Clear polycarbonate student project tanks do not crack or break.

Molded of polycarbonate in either a 6" or 8" depth,

Student project tanks are molded of a shatter-resistant polycarbonate material that resists scratches, dents and clouding-up allowing for long term classroom use. Tanks and lids also naturally resist many acids and oils to maintain a clear student view. Easy to clean and maintain, clear student project tanks will withstand the abuse of daily classroom use.

Available in both 1 and 3 gallon capacities, tanks are supplied with:

- (1) PC clear tank
- (1) PC clear lid
- (1) PC clear drain shelf
- (2) solid Polyamide plugs
- (2) vented Polyamide plugs
- (2) clear plastic / SS clamps

Pop-out plugs are supplied both as solid and vented allowing versatility in use. Included clamps will secure down the lid - keeping all of your critters safe inside!



Polycarbonate can withstand temperatures of -40 to 210°F allowing for flexibility in storage, transport, and use.

Technical Specifications

| 426454 | -0006 | -0008 |
|----------------------|--------|--------|
| Capacity gallon | 1 | 3 |
| Length inches | 10.375 | 12.75 |
| Width inches | 6.375 | 10.375 |
| Depth inches | 6 | 8 |
| Bottom length inches | 8.5 | 11 |
| Bottom width inches | 4.75 | 8.5 |
| Color | Clear | Clear |
| Dishwasher Safe | Yes | Yes |
| Freezer Safe | Yes | Yes |

Ordering Information

| 426454 | Description |
|--------|---|
| -0006 | 1 Gallon Student Project Tank, 6" deep with (4) plugs and (2) clips |
| -0008 | 3 Gallon Student Project Tank, 8" deep with (4) plugs and (2) clips |
| -0000 | Accessory Pack: (4) plugs and (2) clips |

