



"Totally Absorbing" BenchGuard™ Matting



Dynalab®
Labware



Compare and Save!




See product information on reverse of wallet.

BENCHGUARD





Available in sheets and rolls BenchGuard™ Matting



- Highly absorbent paper for bench protection
- Plastic coated underside prevents soak through
- Available in roll or sheet form
- Choice of standard or heavy duty absorption
- BenchGuard Extra has twice the absorption rate at 800 ml/m²

BenchGuard is a proven highly absorbent lab matting that protects lab benches and contains spills. They are both supplied in wallet sheets of 50 each or rolls of 50 meters in a choice of standard or heavy duty absorption. BenchGuard Extra even has twice the absorption rate at 800 ml/m².

504314	Description	Length	Width cm
-0001	Roll dispenser	50 Meters	49
-0002	Extra absorbent roll dispenser	50 Meters	49
-0003	Wallet	60 cm	49
-0004	Extra absorbent wallet	60 cm	49

