



INSTRUCTIONS FOR USE

AZLON® ASPIRATORS & CARBOYS

INTENDED USE

For the storage, transportation and dispensing of larger volumes of media and other liquids.

PREPARATION

- For laboratory use only
- Only the polypropylene product 505634 series is suitable for autoclaving at 121°C – with these items please ensure cap is fully loosened or removed in advance
- Before use with particular liquid, check material resistance and compatibility
- Do not vary intended use of product

HANDLING

- When assembling please ensure stopcock is securely applied and correctly threaded (push, click and tighten)
- Ensure cap is fully engaged to prevent leakage
- When moving follow advised lifting instructions, especially with the larger heavier bottles

STORAGE

- Storage in a dry, clean, well ventilated standard warehouse conditions at temperatures between 15° to 25°C (59°-77°F)

DISPOSAL

- Dispose of in accordance with advised internal company guidelines