



Sample collection and transportation Sterilin Dippas™





Dynalab®
Labware



Sterilin™



Reduce risk of cross-contamination Sterilin Dippas™



- Dippas are designed for sample collection and despatch from one source, thus eliminating the risk of cross contamination
- Dippas are available in transparent polystyrene
- Handles are easily removable, facilitating transport
- All dippas are supplied sterile by irradiation and are individually wrapped to ensure complete sterility
- Blue dippas are highly visible in product samples

507054	Description	Handle Lenth mm	Capacity mL	Cap Material	Height mm	Diameter mm	Case Qty
-0030	Clear	200	30	PP	93.1	30.1	50
-0100	Clear	385	100	Metal	76.4	50.5	100
-0250	Clear	330	250	Metal	120.6	65.1	50
-B030	Blue	200	30	PP	93.1	30.1	50
-B100	Blue	385	100	Metal	76.4	50.5	100
-B250	Blue	330	250	Metal	120.6	65.1	50

