

# Spare Part List

## Automatic Burette

Art. No. 9695-

9695-0001 feed tube, PE, 6 mm Ø


|           |  |       |
|-----------|--|-------|
| 9695-0205 | substitute burette                     | 5 ml  |
| 9695-0210 | substitute burette                     | 10 ml |
| 9695-0215 | substitute burette                     | 15 ml |
| 9695-0225 | substitute burette                     | 25 ml |
| 9695-0250 | substitute burette                     | 50 ml |
| 9695-0216 | substitute burette, without Schellbach |       |

9695-0525 cap without feed tube for 500ml-bottle  
 9695-0528 cap without feed tube for 1000ml-bottle

|           |  |
|-----------|--|
| 9695-0051 | burette holder 500 ml/<br>11 Ø, for 9695-1005, -1010 |
| 9695-0055 | burette hold. 500 ml/15 Ø,<br>for 9695-1015          |
| 9695-0105 | burette hold. 1000 ml/<br>15 Ø, for 9695-1025, -1050 |


9695-0501 narrow-mouth bottle 500 ml without cap  
 0302-1001 narrow-mouth bottle 1000 ml without cap

|           |                     |               |
|-----------|---------------------|---------------|
| 0300-0500 | base for bottle for | 500ml-bottle  |
| 0300-1000 | base for bottle for | 1000ml-bottle |



# Spare Part List Automatic Burette Holder

Art. No. 9695-





# Spare Part List

## Automatic Burette

Art. No. 9695-

| Art. No.  | Article   |  |
|-----------|---|--|
| 9695-0031 | brass pin   |  |
| 9695-0020 | push button                                       |  |
| 9695-0030 | shackle   |  |
| 9695-0032 | micro screw                                       |  |
| 0006-8260 | pressure spring                                   |  |
| 9695-0041 | replacement hose, Latex 100 cm                    |  |
| 9695-0042 | replacement hose, Latex 8 cm                      |  |
| 9695-0040 | exit spout made of polypropylene                  |  |
| 9695-0001 | feed tube, PE, 6 mm Ø                             |  |
| 9695-0525 | cap without feed tube for 500 ml bottle           |  |
| 9695-0528 | cap without feed tube for 1000 ml bottle          |  |
| 9695-0501 | narrow-mouth bottle 500 ml without cap            |  |
| 0302-1001 | narrow-mouth bottle 1000 ml without cap           |  |
| 0300-0500 | base for bottle for 500 ml bottle                 |  |
| 0300-1000 | base for bottle for 1000 ml bottle                |  |
| 9695-0205 | substitute burette glass 5 ml                     |  |
| 9695-0210 | substitute burette glass 10 ml                    |  |
| 9695-0215 | substitute burette glass 15 ml                    |  |
| 9695-0225 | substitute burette glass 25 ml                    |  |
| 9695-0250 | substitute burette glass 50 ml                    |  |
| 9695-0051 | burette holder 500 ml/11 Ø, for 9695-1005, -1010  |  |
| 9695-0055 | burette holder 500 ml/15 Ø, for 9695-1015         |  |
| 9695-0105 | burette holder 1000 ml/15 Ø, for 9695-1025, -1050 |  |