

Dispenser Bottles

with Adjustable Volume

Kartell

Dispenser Bottle with Adjustable Volume LDPE/PMP

- LDPE bottle with PMP dispenser

The graduated clear PMP measuring cup slides up and down on the vertical polyethylene tube for easy adjustment. The variation in height determines the volume to be dispensed. The cup is filled by squeezing the bottle and forcing the liquid up the tube. When the bottle pressure is released, the excess will be drawn back into the bottle and only the desired volume will remain.



	206305	206315	206325
Bottle Size ml	250	500	1000
Bottle Diameter mm	60	75	94
Cup Size ml	5-25	5-25	5-25
Height mm	282	324	354
Case Quantity	10	10	10

	206335	206345	206355
Bottle Size ml	250	500	1000
Bottle Diameter mm	60	75	94
Cup Size ml	5-50	5-50	5-50
Height mm	298	340	370
Case Quantity	10	10	10