



# Laboratory Organizer Spill Stand Carboy Spill Containment

**Dynalon**  
Labware  
quality . value . service



Designed for laboratory carboys and other large containers with spigots to contain liquid spills.

**Custom sizes  
available  
upon request.**



175 Humboldt Street, Suite 300  
Rochester, NY 14610  
Telephone 800-334-7585  
Fax 585-334-0241  
[www.dynalon.com](http://www.dynalon.com)



# Dynalon Labware is your source for safety products!

Built-in lower bin for convenient storage of labware or small accessories.

- Trays are fabricated with 1/4" thick chemically resistant HDPE. All welded construction featuring a 1/2" skid resistant top shelf
- Provides a safer way to dispense chemicals and other liquids in the lab
- Protects cabinets, laboratory benches and floors from liquid spills
- Each tray contains two convenient drain ports and 1/2" tubulation connection on lower trough
- Built-in side handles for convenient transport
- 8" shelf height to accommodate up to a 5000mL beaker
- All trays are fabricated and tested in our Rochester, New York facility
- Custom sizes available upon request.



#### 107364 Carboy Spill Tray Specifications:

Height: 13" / 330.2 mm

Width: 24.5" / 622.3 mm

Depth: 25" / 635 mm

Shelf: 24"W x 12"D / 609 x 305 mm

Tray (lower trough) capacity: 2.5 Gal / 10 L

Weight: 16 lbs / 7.25 Kg

**Dynalon®**  
Labware

quality . value . service

175 Humboldt Street, Suite 300

Rochester, NY 14610

Telephone 800-334-7585

Fax 585-334-0241

[www.dynalon.com](http://www.dynalon.com)