

# *Dippas<sup>®</sup> sterile sampler* *and test container all in one!*




Dippas are individually wrapped to ensure complete sterility.



Available in either transparent polystyrene or shatterproof polypropylene in three different capacities.

Blue Dippas are highly visible in product samples.



Designed for sample collection and dispatch from one source, thus eliminating the risk of cross contamination.



Closures are supplied in either PP or Metal in PS Dippas and in PE or PP in Polypropylene dippers.



Handles are easily removed, enabling the sample to be sent away for analysis.



Sample collection and transportation to the laboratory in a single container.