


## DispoLadle

DispoLadle is a sampling scoop specially designed for disposable use. The scaled scoop with embossing can be used for measuring. It shows contents in increments of 10 ml, up to 200 ml. An ergonomic handle and a spout on both sides for right or left-hand operation make handling simple and easy.

- ▶ Production, assembly & packaging in Class 7 clean room (10,000)
- ▶ 200 ml scoop - vertical
- ▶ External Ø 92 mm (3.62")
- ▶ Ergonomic handle
- ▶ HDPE - food grade
- ▶ Clean room manufactured
- ▶ Individually packaged
- ▶ Optionally sterilized by gamma rays



### DispoLadle

Length mm	Model	Pieces per pack	Item no.		
500	nonsterile	20	853933-6640	N4	
500	sterile	20	853933-6641	N4	
1000	nonsterile	20	853933-6650	N4	
1000	sterile	20	853933-6651	N4	

