

Flip & Spray™ Bottle PE/PP

The spray bottle that even works **upside down!**

These unique 250, 500 and 1000mL reusable polyethylene Spray Bottles with a polypropylene sprayer have excellent chemical resistance and will perform upright or inverted. The limber silicone tubing attached to a stainless steel weight with a PP filter remains immersed in the liquid in the upside down position and allows the product to work. This unique feature lets a user spray conveniently from any direction while reaching many inaccessible areas. Adjust the spray knob to obtain a fine mist to a steady stream.

- Uninterrupted spraying
- Parts in contact with media made of PP, PE, silicon and stainless steel
- Nozzle diameter 0.6 mm and continuously adjustable
- Amount sprayed per stroke: 1.2 ml ± 0.1

Product Specifications

805104	-0250	-0500	-1000
Bottle Capacity mL	250	500	1000
Height w Closure mm	220	240	295
Neck ID mm	18	18	17.5



FREE

Buy Any 3
Flip n Spray Bottles,
receive a
**FREE write-on
500mL Wash Bottle!**

Please list promo code **FSB18** on purchase order.

Dynalon®
Labware
quality . value . service