

Mini Solvent Pump

Reliable filling of small quantities



- Ideal for 1 to 5 Liter containers
- Connecting thread GL45
- Parts in contact with medium: PTFE and stainless steel
- PTFE tubing (6x8mm) 0.4m
- Includes grounding cable 2m
- Screw thread adapter available as accessory for S 40, GL 38 and GL 32 threads



Dynalon
Labware
quality . value . service

175 Humboldt Street, Suite 300
Rochester, NY 14610
Tel 800-334-7585
Fax 585-334-0241
CS@dyna-labare.com

www.dynalon.com

Free Assortment of 500mL GHS Wash Bottles with the purchase of each pump!



NEW Mini Solvent Pump SS/PTFE GL45

Precise dispensing of liquids.

The Mini Solvent Pump is a high-quality hand pump for use with liquids and flammable media such as acetone, ethanol, isopropanol and methanol. Use with small containers such as bottles, canisters, carboys and tanks with GL45 threads, from 1 liter up to ~10 liters in size. Accessory thread adapters are available that allow this pump to fit on S40, GL38 and GL32 threads.

The pump works according to the excess pressure principle. Through use of the pump ball and the PTFE tubing, liquid is forced from the container due to the pressure build up. The pump includes a pressure relief valve for safety purposes.

Parts that come into contact with the media are stainless steel (1.4301 and 1.4404) and PTFE tubing cut to size; both have excellent chemical resistance. Also included is a grounding cable that easily inserts into the pump housing, for the purpose of diverting any electrostatic charges.

The Mini Solvent Pump is compact and very well balanced. Dosing of liquids occurs through the valve with an auto-stop function. It is very precise with no dripping even into small containers such as test tubes. Tubing can even be attached to the discharge spout to fill containers at a further distance. The high quality materials ensure neat and reliable filling throughout a very long service life.

Mini Solvent Pump includes:

- PTFE tubing 0.4 meters
- Grounding cable, 2 meters

#856033-4000

Mini Solvent Pump SS/PTFE GL45 \$695.00

Accessories:

#0006-4000	Thread Adapter PTFE GL45-S40	\$53.58
#0006-4002	Thread Adapter ETFE GL45-GL38	\$53.58
#0006-4003	Thread Adapter ETFE GL45-GL32	\$53.58
#0006-9506	Tubing Cutter	\$13.78
#8804-0608	PTFE Tubing, 6 x 8mm, 5 meter Roll	\$74.09



Receive a Free Assortment of
500mL GHS Wash Bottles
with the purchase of each pump!
\$40 value

Good on orders received through August 31, 2017.