


2017 NEW Plastic Lab Supplies



NEW Test Tube Racks, PP

Kartell

Lightweight, compact Polypropylene test tube rack allows contents of tube to be seen from top to bottom. Specifically designed test tube rack for titration work and extremely useful for viewing colors and precipitates. May be immersed in water bath and can withstand temperatures up to 120C.

	259315	259325	259335
Places	8	9	10
For Tube Diameter mm	20	18	16
Length inches	8.66	8.66	8.66
Width inches	2.12	2.17	2.17
Height inches	2.58	2.56	2.46
Price Each \$	10.61	10.50	10.70

NEW Wide Orifice Pipette, PP

- Sampling small volumes of viscous substances
- 10mL capacity with 0.1mL sub divisions

Translucent polypropylene. Resistant to breakage. Blue graduations, ring mark and inscriptions. Use with pipettor. Recommended cleaning with mild alkaline detergent up to 60°C.

Note: Exposure to temperatures above 60°C may cause volume changes. Additional sizes available upon request.

102774-0010 \$15.06 each



NEW Graduated Bottles, PC with PP Closures

- Autoclavable and leakproof
- Glass clear, durable and highly break resistant
- Graduated in metric and English;
- Very high quality

Boston Round bottles, Narrow Mouth. Use with biologicals, liquids, reagents, media, buffers and aqueous solutions. Manufactured in Class 10K cleanroom. Available in packs or shrink wrap trays.

	-0300	-0500	-1000
601524			
Size mL	300	500	1000
Closure mm	38-430	38-430	38-430
Pack Quantity	6	6	4
Pack Price \$	\$45.30	\$56.82	\$44.76
Case Quantity	24	12	12
Case Price \$	\$181.20	\$113.64	\$134.20

175 Humboldt Street, Suite 300
Rochester, NY 14610

Tel 800-334-7585
Fax 585-334-0241
www.dynalon.com

NEW Flip & Spray™ Bottle PE/PP

- Uninterrupted spraying
- Parts in contact with media made of PP, PE, silicon and stainless steel
- Nozzle diameter 0.6 mm and continuously adjustable
- Amount sprayed per stroke: 1.2 ml ± 0.1

These unique 250, 500 and 1000mL reusable polyethylene Spray Bottles with a polypropylene sprayer have excellent chemical resistance and will perform upright or inverted. The limber silicone tubing attached to a stainless steel weight with a PP filter remains immersed in the liquid in the upside down position and allows the product to work. This unique feature lets a user spray conveniently from any direction while reaching many inaccessible areas. Adjust the spray knob to obtain a fine mist to a steady stream.

	-0250	-0500	-1000
805104			
Bottle Capacity mL	250	500	1000
Height w Closure mm	220	240	295
Neck ID mm	18	18	17.5
Price Each \$	16.03	17.50	18.42

NEW Carboy Spill Containment Tray and Lab Organizer, HDPE

- Designed for laboratory carboys or other large containers with spigots
- Trays are fabricated with 1/4" thick chemically resistant HDPE.
- All welded construction featuring a 1/2" skid resistant top shelf.
- Provides a safer way to dispense chemicals or other liquids in the lab
- Protects cabinets, laboratory benches and floors from liquid spills
- Each tray contains two convenient drain ports and 1/2" tubulation connection on lower trough
- Built-in lower bin for convenient storage of labware or accessories
- Built-in side handles for convenient transport
- 8" shelf height to accommodate up to a 5000mL beaker
- All trays are fabricated and tested in our Rochester, New York facility

Height: 13" / 330.2 mm
Width: 24.5" / 622.3 mm
Depth: 25" / 635 mm
Shelf: 24"W x 12"D / 609 x 305 mm
Tray (lower trough) capacity: 2 1/2 Gal / 10 L
Weight: 16 lbs / 7.25 Kg

107364 Carboy Spill Tray \$357.54 each



NEW GHS Wash Bottle, LDPE, Bleach

- 500mL Capacity
- Multi-lingual wide neck wash bottle
- Opaque white LDPE for light sensitive sodium hypochlorite / bleach
- Indelibly printed in three languages - English, French and Spanish
- Features the new Globally Harmonized Symbols, NFPA diamond, risk phrases and safety information
- Fitted with color coded one-piece yellow PP closure
- Closure features precision jet-tip aperture for fine control
- Dimensions 163 h x 75mm diameter



506975 \$36.20 pack of 5



NEW Dynalon® Beakers with Handle in Metric and English, PP

- Made from translucent high clarity polypropylene - food grade resin
- Graduated in Metric and English in attractive blue print with scales on both sides - easy to read for right and left handed use
- Tapered and space saving; nest for storage
- Durable handle and thumb grip for easy dispensing
- Handle allows for drainage when inverted
- Maximum operating temperature 125°C

	-0500	-1000
326525		
Capacity ml	500	1000
Sub Div. ml	10	10
Sub Div. oz	1	1
Diameter mm	97	127
Height mm	138	168
Price Each \$	10.20	13.25

NEW Stopcock Tubing 2-Way

- PP/PE: for industrial use, red/white
- Operating temperature: PP/PE +5 up to +40 °C
- Resistant to compression max. 1 bar or 14.5 psi
- With mounting plate for easy installation

	-0060	-0080	-0100	-0120
886054				
Material	PP / PE	PP / PE	PP / PE	PP / PE
For diameter mm	5 - 7	7 - 9	9 - 11	11 - 13
Nominal width mm	4	6	8	8
Price Each \$	15.65	18.04	20.69	20.69



The spray bottle that even works upside down!



NEW Straight-Side Containers PETE with White PP Closures

Transparent Straight-Side Containers of PETE (Polyethylene Terephthalate) with white polypropylene liner less screw closures. Store pastes, samples, powders, specimens etc. in industrial and lab applications. PETE is durable, shatterproof and has excellent gas barrier properties with H₂O, O₂ and CO₂. These jars may be used below freezing to a maximum of 70°C.

426315	-0004	-0008	-0016	-0032
Capacity mL	125	250	500	1000
Diameter mm	73	73	93.5	92.4
Height with closure mm	46.4	79.5	86	173
Pack Quantity	12	6	6	6
Pack Price \$	28.80	18.30	25.92	39.60

NEW Petri Dishes, PS

- Transparent Polystyrene and optically clear
- Single use and stackable
- Sterile by irradiation (-0000 is gamma irradiated; -0004 is E-Beam irradiated)
- Consistent flatness allows for even media distribution
- Certificate of Sterilization upon request
- Traceable by batch number

402024	-0000	-0004
Size mm (nominal)	100	100
Color	Clear	Clear
Vents	None	Four
Dishes per Sleeve	20	25
Quantity per Box	480	500
Case Price	\$121.44	\$126.48



NEW Autoclavable Deep Tote, PP

Tough 12 liter autoclavable polypropylene tray for a variety of uses. Can hold 18 x 500 ml or 12 x 1 liter bottles. Tapered and stackable to save storage space. Product Dimensions (l x w x h): 20 x 13.8 x 4.1"

210505 \$31.50 each



NEW Bucket with Graduations and Spout, LDPE

These LDPE buckets with easy pour spout are graduated internally in liters (TC-to contain). The steel bale handle has a comfortable grip. Suitable for food and beverage products with good to excellent chemical resistance especially at room temperature. Use up to 80°C or well below freezing.

213185	-0009	-0012	-0017
Capacity liters	9	12	17
Capacity gallons	2.37	3.17	4.49
O.D. inches	9.84	11.02	12.2
Height inches	11.02	11.81	13.38
Price Each \$	16.75	22.15	32.35



NEW Wash Bottle, Autoclavable, PP

- Autoclavable at 121°C
- PP 500mL bottle fitted with one-piece PP closure (34mm neck ID)
- Closures feature precision jet-tip aperture for fine control
- Dimensions 163 height x 75mm diameter

506985 \$44.95 pack of 5



NEW Storage and Transport Trays, PVC

Useful for Lab and Desk Items. Very sturdy. Use for carrying items such as stir bars, connectors etc. Useful for knick knacks - rulers, erasers, pencils etc, in desk drawers. Central compartment of five cavity tray can be used for pipettes.

209304	-0005	-0012
Cavity	5	12
Cavity dimensions mm	90 x 90	100 x 185 (4) 100 x 385 (1)
Length mm	303	304
Width mm	403	404
Height mm	63	64
Price Each \$	53.85	58.25



NEW 30 Liter Bio-bin® Spill Containment Trays, PP

- Ideal spill trays for autoclaving waste disposal containers such as the Dynalon Bio-bins
- Autoclave tray for other plastic or glass laboratory items
- Help contain accidental chemical spills in the laboratory
- Material: Polypropylene (PP)
- 3/16" double-welded construction
- Tested in the USA



797303	-6310	-6320	-6330
Description	Upright	Side Flaps Closed	Side Flaps Open
Volume gallons	2	3	4
Length inches	12	12	12
Width inches	12	18	26
Height inches	3	3	3
Price Each \$	93.83	115.80	145.11

Dry Heat Sterilization Tape, 50 Meters



Use these tapes to enhance confidence that sterilization has occurred. Color changes from white to brown. These tapes conform with ISO 1140/2005.

202075 \$52.85 per roll

NEW Tuff Vials, PC

Highly break resistant tough polycarbonate vials with polypropylene leakproof closures. Glass clear allows viewing of valuable contents, liquid or solid. Reusable and autoclavable. 180 mL can also be used as an outer protective container. Sold as each, Packs of 6 or Tray of 12.



5mL and 10mL are also ideal for leakproof shipping of reagents and buffer solutions, for packaging of serum, veterinary, pharmaceuticals, cosmetic and nutraceutical extracts. Use with most aromatic compounds and essential oils. Available in Packs of 20 or shrink wrap trays of 100.

626434	-0005	-0010	-0180
Size mL	5	10	180
Closure mm	20-415	20-415	54
Pack Quantity	20	20	6
Pack Price \$	\$16.00	\$17.60	\$30.96
Case Quantity	100	100	12
Case Price \$	\$80.00	\$88.00	\$63.00



NEW Tray, Sample Bottle, PVC

A strong PVC tray for holding sample bottles with base diameter not greater than 22 mm. Ideal for 7 ml (1/4oz) Bijou bottles or specimen tubes. Seen in photo are new blue closed Tuff vials (below), 5 & 10mL, # 626434. Each of nine rows can hold 12 to 15 of such bottles making a total capacity of 108 or more. Dimensions: 350 x 300 x 40mm.

259214 \$51.20 each



NEW Tray, Pipette 4 Place, PVC

Designed to hold ~30 pipettes. 4 compartments accommodate pipettes of 1-2-5-10 ml capacities. Dimensions: 510 x 350 x 105mm.

211004 \$35.96 each



NEW Stackable Tanks, HDPE

Can be stacked to minimize storage space. High chemical resistance to organic and inorganic acids. Robust with lifting lugs and strengthening ribs to ensure long life.

219274	-0020	-0046
Capacity liters	20	46
Length inches	16.3	23.4
Width inches	12.4	14.4
Height inches	7.9	10
Price Each \$	45.70	85.30

