

Titrating Burette Automatic

Available in 5 sizes from 5mL to 50mL

This automatic titrating burette is a universal aid in the laboratory, electroplating industry, water treatment plants, for field analyses and numerous other applications. Available in 5 sizes from 5mL to 50mL.

- Measured "to deliver"
- Burette made of Schott graduated borosilicate glass with blue ring graduations and Schellbach stripe for easy to read measurements
- Graduations in accordance with Category B (Tolerances according to DIN ISO 385)
- Exact dispensing of precise discharge jet or drop via the micro screw
- Left or right hand operation for ease of use
- Sturdy stand for support
- 1 liter low density polyethylene (LDPE) chemically resistant bottle

Product Specifications

896954	-1005	-1010
Burette Capacity mL	5	10
Bottle Capacity mL	500	500
Division mL	0.05	0.05
Precision mL	± 0.015	± 0.030
Height cm	48	57
Price Each \$	118.49	118.49

896954	-1015	-1025
Burette Capacity mL	5	25
Bottle Capacity mL	500	1000
Division mL	0.10	0.10
Precision mL	± 0.045	± 0.045
Height cm	49	61
Price Each \$	118.49	118.49

896954	-1050
Burette Capacity mL	50
Bottle Capacity mL	1000
Division mL	0.10
Precision mL	± 0.045
Height cm	87
Price Each \$	138.42



Dynalon®
Labware
quality . value . service

175 Humboldt Street, Suite 300
Rochester, NY 14610
Tel 800-334-7585
Fax 585-334-0241
CS@dyna-labare.com