



Dynalab Corp.

Dynalon Labware

175 Humboldt Street, Suite 300
Rochester, New York 14610 USA

Tel 585.334.2060 800.828.6595
Fax 585.334.0241
www.DynalabCorp.com

Certificate of Compliance

"EQUAL Y" Tubing Connector, Polypropylene

Dynalon Products Covered:

227005, 227025, 227045, 227065, 227085 and 227105

Regarding our "EQUAL Y" Tubing Connectors, I wish to inform you that these are manufactured in Italy from autoclavable Polypropylene resin. Batch information such as expiration date is included on Kartell package labels when full cases are ordered.

"EQUAL Y" Tubing Connectors have excellent chemical resistance and withstands flows to 20 p.s.i. Serrated connecting surfaces to prevent tubing slipping. Polypropylene "EQUAL Y" Tubing Connectors will withstand temperatures up to 135°C.

Compliance Information (provided by manufacturer):					
BSE/TSE Free	Yes	FDA Compliant	Yes	RoHS 3	Yes
California Proposition 65	Yes	Reach	Yes		
Mfg. Part Numbers	527, 528, 529, 530, 531, 532				
Mfg. Certifications	UNI EN ISO 9001:2015		UNI EN ISO 14001:2015		

Please advise if you require any further information.

Regards,

Robert Pfeil
Vice President Operations / Quality Manager
Dynalab Corp | Dynalon Labware