



**Dynalab Corp.**

**Dynalon Labware**

175 Humboldt Street, Suite 300  
Rochester, New York 14610 USA

Tel 585.334.2060 800.828.6595  
Fax 585.334.0241  
www.DynalabCorp.com

## Certificate of Compliance

### "L" Tubing Connector, Polypropylene

Dynalon Products Covered:

**227165, 227175, 227185, 227195, 227205, 227215 and 227225**

Regarding our "L" Tubing Connectors, I wish to inform you that these are manufactured in Italy from autoclavable Polypropylene resin. Batch information such as expiration date is included on Kartell package labels when full cases are ordered.

"L" Tubing Connectors have excellent chemical resistance and withstands flows to 20 p.s.i. Serrated connecting surfaces to prevent tubing slipping. Polypropylene "L" Tubing Connectors will withstand temperatures up to 135°C.

Compliance Information (provided by manufacturer):					
BSE/TSE Free	Yes	FDA Compliant	Yes	RoHS 3	Yes
California Proposition 65	Yes	Reach	Yes		
Mfg. Part Numbers	1250, 1251, 1252, 1253, 1254, 1255, 1256				
Mfg. Certifications	UNI EN ISO 9001:2015		UNI EN ISO 14001:2015		

Please advise if you require any further information.

Regards,

Robert Pfeil  
Vice President Operations / Quality Manager  
Dynalab Corp | Dynalon Labware