



**Dynalab Corp.**

**Dynalon Labware**

175 Humboldt Street, Suite 300  
Rochester, New York 14610 USA

Tel 585.334.2060 800.828.6595  
Fax 585.334.0241  
www.DynalabCorp.com

## Certificate of Compliance

### Kartell Connector Stopper 29/32, Polypropylene

Dynalon Product Covered:

**228205**

Regarding our Polypropylene Connector Stopper 29/32, I wish to inform you that these are manufactured in Italy with an autoclavable Polypropylene resin.

Kartell plastic tubing connector stoppers with tubing connections. Standard taper 1:10 ensuring perfect sealing action. Hose adapters are Ø mm 8 (~5/16") - Length 20 mm.

Compliance Information (provided by manufacturer):					
BSE/TSE Free	Yes	FDA Compliant	Yes	RoHS 3	Yes
California Proposition 65	Yes	Reach	Yes		
Mfg. Part Number	1816				
Mfg. Certifications	UNI EN ISO 9001:2015		UNI EN ISO 14001:2015		

Please advise if you require any further information.

Regards,

Robert Pfeil  
Vice President Operations / Quality Manager  
Dynalab Corp | Dynalon Labware