


Technical Datasheet

Article - No:	TDS-WMPB015H	Revision: A0																							
Remarks	single use only																								
Specification																									
Product name	WMPB015H,15ml Packaging Bottle,HDPE																								
Production conditions	Class 100,000 clean room production																								
	No release agent is used																								
Classification	CE/IVDD																								
Dimensions	<div></div> <table><tr><th>location name</th><th>material</th><th>O. D. D1 (mm)</th><th>Height H1 (mm)</th><th>I. D. D2 (mm)</th><th>Cap O.D D3 (mm)</th><th>Cap H H2 (mm)</th><th>Total H (mm)</th></tr><tr><td rowspan="2">15ml</td><td>PP</td><td>φ24.8±0.2</td><td>56.0±0.2</td><td>φ13.8±0.2</td><td>φ24.2±0.2</td><td>14.5±0.2</td><td>57.2±0.5</td></tr><tr><td>PE</td><td>φ24.8±0.2</td><td>56.1±0.2</td><td>φ13.8±0.2</td><td>φ24.2±0.2</td><td>14.5±0.2</td><td>57.3±0.5</td></tr></table>		location name	material	O. D. D1 (mm)	Height H1 (mm)	I. D. D2 (mm)	Cap O.D D3 (mm)	Cap H H2 (mm)	Total H (mm)	15ml	PP	φ24.8±0.2	56.0±0.2	φ13.8±0.2	φ24.2±0.2	14.5±0.2	57.2±0.5	PE	φ24.8±0.2	56.1±0.2	φ13.8±0.2	φ24.2±0.2	14.5±0.2	57.3±0.5
	location name	material	O. D. D1 (mm)	Height H1 (mm)	I. D. D2 (mm)	Cap O.D D3 (mm)	Cap H H2 (mm)	Total H (mm)																	
15ml	PP	φ24.8±0.2	56.0±0.2	φ13.8±0.2	φ24.2±0.2	14.5±0.2	57.2±0.5																		
	PE	φ24.8±0.2	56.1±0.2	φ13.8±0.2	φ24.2±0.2	14.5±0.2	57.3±0.5																		
Colour	Clear																								
Nominal Capacity	15ml																								
Material Characteristics	bottle:High Density Polyethylene(HDPE)																								
	cap:polypropylene (PP)																								
Performance parameters	Negative pressure:-70kpa~-90kpa																								
	Ethanol loss rate:≤0.5%																								
	water loss rate:≤0.2%																								
Purity Grade	free of RNase and DNase																								
	free of pyrogen																								
	free of PCR inhibitors																								
	free of human DNA																								
	free of endotoxins																								
Shelf - Life	36 months at controlled storage conditions																								
Packaging	pcs per case: 1,500																								
Storage	dry condition:40%-60% rel. humidity																								
Approved by quality department																									