


## Technical Datasheet

Article - No:	TDS-WMPB060H	Revision: A0																							
Remarks	single use only																								
Specification																									
Product name	WMPB060H,60ml Packaging Bottle,HDPE																								
Production conditions	Class 100,000 clean room production																								
	No release agent is used																								
Classification	CE/IVDD																								
Dimensions	<div></div> <table><tr><th>location name</th><th>material</th><th>O. D. D1 (mm)</th><th>Height H1 (mm)</th><th>I. D. D2 (mm)</th><th>Cap O.D D3 (mm)</th><th>Cap H H2 (mm)</th><th>Total H (mm)</th></tr><tr><td rowspan="2">60ml</td><td>PP</td><td>φ39.0±0.3</td><td>81.5±0.3</td><td>φ21.1±0.3</td><td>φ32.2±0.2</td><td>15.5±0.2</td><td>83.5±0.5</td></tr><tr><td>PE</td><td>φ38.6±0.3</td><td>81.65±0.3</td><td>φ20.9±0.3</td><td>φ32.2±0.2</td><td>15.5±0.2</td><td>83.65±0.5</td></tr></table>		location name	material	O. D. D1 (mm)	Height H1 (mm)	I. D. D2 (mm)	Cap O.D D3 (mm)	Cap H H2 (mm)	Total H (mm)	60ml	PP	φ39.0±0.3	81.5±0.3	φ21.1±0.3	φ32.2±0.2	15.5±0.2	83.5±0.5	PE	φ38.6±0.3	81.65±0.3	φ20.9±0.3	φ32.2±0.2	15.5±0.2	83.65±0.5
	location name	material	O. D. D1 (mm)	Height H1 (mm)	I. D. D2 (mm)	Cap O.D D3 (mm)	Cap H H2 (mm)	Total H (mm)																	
60ml	PP	φ39.0±0.3	81.5±0.3	φ21.1±0.3	φ32.2±0.2	15.5±0.2	83.5±0.5																		
	PE	φ38.6±0.3	81.65±0.3	φ20.9±0.3	φ32.2±0.2	15.5±0.2	83.65±0.5																		
Colour	Clear																								
Nominal Capacity	60ml																								
Material Characteristics	bottle:High Density Polyethylene(HDPE)																								
	cap:polypropylene (PP)																								
Performance parameters	Negative pressure:-70kpa~-90kpa																								
	Ethanol loss rate:≤0.5%																								
	water loss rate:≤0.2%																								
Purity Grade	free of RNase and DNase																								
	free of pyrogen																								
	free of PCR inhibitors																								
	free of human DNA																								
	free of endotoxins																								
Shelf - Life	36 months at controlled storage conditions																								
Packaging	pcs per case:800																								
Storage	dry condition:40%-60% rel. humidity																								
Approved by quality department																									