


Technical Datasheet

Article - No:	TDS-WMPB125	Revision: A0																							
Remarks	single use only																								
Specification																									
Product name	WMPB125,125ml Packaging Bottle,HDPE																								
Production conditions	Class 100,000 clean room production																								
	No release agent is used																								
Classification	CE/IVDD																								
Dimensions	<div></div> <table><tr><th>location name</th><th>material</th><th>O. D. D1 (mm)</th><th>Height H1 (mm)</th><th>I. D. D2 (mm)</th><th>Cap O.D D3 (mm)</th><th>Cap H H2 (mm)</th><th>Total H (mm)</th></tr><tr><td rowspan="2">125ml</td><td>PP</td><td>φ50.7±0.4</td><td>95.7±0.4</td><td>φ28.0±0.4</td><td>φ42.2±0.3</td><td>15.5±0.3</td><td>97.8±0.8</td></tr><tr><td>PE</td><td>φ50.8±0.4</td><td>96.2±0.4</td><td>φ27.8±0.4</td><td>φ42.2±0.3</td><td>15.5±0.3</td><td>98.3±0.8</td></tr></table>		location name	material	O. D. D1 (mm)	Height H1 (mm)	I. D. D2 (mm)	Cap O.D D3 (mm)	Cap H H2 (mm)	Total H (mm)	125ml	PP	φ50.7±0.4	95.7±0.4	φ28.0±0.4	φ42.2±0.3	15.5±0.3	97.8±0.8	PE	φ50.8±0.4	96.2±0.4	φ27.8±0.4	φ42.2±0.3	15.5±0.3	98.3±0.8
location name	material	O. D. D1 (mm)	Height H1 (mm)	I. D. D2 (mm)	Cap O.D D3 (mm)	Cap H H2 (mm)	Total H (mm)																		
125ml	PP	φ50.7±0.4	95.7±0.4	φ28.0±0.4	φ42.2±0.3	15.5±0.3	97.8±0.8																		
	PE	φ50.8±0.4	96.2±0.4	φ27.8±0.4	φ42.2±0.3	15.5±0.3	98.3±0.8																		
Colour	Clear																								
Nominal Capacity	125ml																								
Material Characteristics	bottle:High Density Polyethylene(HDPE)																								
	cap:polypropylene (PP)																								
Performance parameters	Negative pressure:-70kpa~-90kpa																								
	Ethanol loss rate:≤0.5%																								
	water loss rate:≤0.2%																								
Purity Grade	free of RNase and DNase																								
	free of pyrogen																								
	free of PCR inhibitors																								
	free of human DNA																								
	free of endotoxins																								
Shelf - Life	36 months at controlled storage conditions																								
Packaging	pcs per case:500																								
Storage	dry condition:40%-60% rel. humidity																								
Approved by quality department																									