


Technical Datasheet

Article - No:	TDS-WMPB250	Revision: A0																							
Remarks	single use only																								
Specification																									
Product name	WMPB250,250ml Packaging Bottle,HDPE																								
Production conditions	Class 100,000 clean room production																								
	No release agent is used																								
Classification	CE/IVDD																								
Dimensions	<div></div> <table><tr><th>location name</th><th>material</th><th>O. D. D1 (mm)</th><th>Height H1 (mm)</th><th>I. D. D2 (mm)</th><th>Cap O.D D3 (mm)</th><th>Cap H H2 (mm)</th><th>Total H (mm)</th></tr><tr><td rowspan="2">250ml</td><td>pp</td><td>φ60.9±0.4</td><td>127.0±0.4</td><td>φ33.1±0.4</td><td>φ46.4±0.4</td><td>16.8±0.4</td><td>129.8±0.8</td></tr><tr><td>PE</td><td>φ61.7±0.4</td><td>127.6±0.4</td><td>φ33.0±0.4</td><td>φ46.4±0.4</td><td>16.8±0.4</td><td>130.8±0.8</td></tr></table>		location name	material	O. D. D1 (mm)	Height H1 (mm)	I. D. D2 (mm)	Cap O.D D3 (mm)	Cap H H2 (mm)	Total H (mm)	250ml	pp	φ60.9±0.4	127.0±0.4	φ33.1±0.4	φ46.4±0.4	16.8±0.4	129.8±0.8	PE	φ61.7±0.4	127.6±0.4	φ33.0±0.4	φ46.4±0.4	16.8±0.4	130.8±0.8
	location name	material	O. D. D1 (mm)	Height H1 (mm)	I. D. D2 (mm)	Cap O.D D3 (mm)	Cap H H2 (mm)	Total H (mm)																	
250ml	pp	φ60.9±0.4	127.0±0.4	φ33.1±0.4	φ46.4±0.4	16.8±0.4	129.8±0.8																		
	PE	φ61.7±0.4	127.6±0.4	φ33.0±0.4	φ46.4±0.4	16.8±0.4	130.8±0.8																		
Colour	Clear																								
Nominal Capacity	250ml																								
Material Characteristics	bottle:High Density Polyethylene(HDPE)																								
	cap:polypropylene (PP)																								
Performance parameters	Negative pressure:-70kpa~-90kpa																								
	Ethanol loss rate:≤0.5%																								
	water loss rate:≤0.2%																								
Purity Grade	free of RNase and DNase																								
	free of pyrogen																								
	free of PCR inhibitors																								
	free of human DNA																								
	free of endotoxins																								
Shelf - Life	36 months at controlled storage conditions																								
Packaging	pcs per case:250																								
Storage	dry condition:40%-60% rel. humidity																								
Approved by quality department																									