


Good for ~5,000 pipettings MiniFIX Micropipettes



Dynalab®
Labware

• AUTOCLAVABLE

- Cost effective
- Fixed volume from 5uL to 1000uL
- Ideal in schools for routine applications where speed and precision is required
- Color-coded button for easy volume identification
- Accommodates all common tips



MiniFIX line of fixed volume mini micropipettes consists of ten pipettes covering the range from 5 uL to 1000 uL.

These pipettes feature a compact design, with a length of approximately 13cm, and weight of less than 10g. They are single blister packed, with calibration and performance certification. Each MiniFIX is marked with its volume, and the push button is color coded for easy "volume" identification.



240234	Vol. uL	Acc. %	C.V. %
-0005	5	±1.5	±1.0
-0010	10	±1.0	±1.0
-0020	20	±0.5	±0.5
-0025	25	±0.5	±0.5
-0050	50	±0.4	±0.3
-0100	100	±0.3	±0.3
-0200	200	±0.3	±0.2
-0250	250	±0.3	±0.2
-0500	500	±0.3	±0.2
-1000	1000	±0.3	±0.2



MiniFIX Rack

Made of transparent acrylic, the MiniFIX Microliter Pipette Rack will hold up to 5 microliter pipettes. Helps to organize your lab and keep the pipettes at your finger tips.

#140234