No Vent Petri Dish



Dynalon® Labware

quality . value . service

175 Humboldt Street, Suite 300 Rochester, NY 14610 Tel 800-334-7585 Fax 585-334-0241

NEW #402024-0000

- No vent 92 x 16mm stackable petri dish
- Transparent top and bottom
- Disposable (single use) polystyrene petri dis
- Optically clear
- Gamma sterilized
- Contents sterile if sleeve is unopened
- Consistent flatness allows for even media distribution
- Certificates of Sterilization Processing upon request
- Traceable by batch number

Dishes/	Dishes/	Dishes/	Box Dimensions			Box	Box
Sleeve	Box	Pallet	L W H			Cubic	Weight
20	480	9600	23 584	15.24" 387mm	14.17 360	2.88ft³ 0.0814m³	24.62lbs 9.19 kg


Stackable Petri Dish



Dynalon® Labware

quality . value . service

175 Humboldt Street, Suite 300 Rochester, NY 14610 Tel 800-334-7585 Fax 585-334-0241



Dishes/	Dishes/	Boxes/	Box Dimensions			Box	Box
Sleeve	Box	Pallet	L W H			Cubic	Weight
25	500	30	19.25 48.9	15.5" 39.37cm	16.0 40.64	2.76ft³ 0.08m³	20.86lbs 9.46 kg

NEW #402024-0004

- Standard 100 x 15mm stackable petri dish
- Clear top and bottom viewing surfaces
- Disposable (single use) polystyrene petri dish for general microbiological use
- Molded from medical-grade virgin polystyrene
- · Optically clear
- Radiation sterilized
- Contents sterile if sleeve is unopened
- 4 Ventilation ribs/vents allow for free air circulation and reduce condensation during incubation
- Consistent flatness allows for even media distribution
- Flared lid skirt and squared corners for easy one-hand operation, especially with gloves
- Certificates of Sterilization Processing upon request
- Traceable by lot and SKU number
- Certificates of Compliance available available upon request



Lid designed for optimum stack stability



Manufactured from virgin-grade polystyrene for optimum clarity and quality



Ventilation ribs allow for free air circulation and reduce condensation during incubation