




**NEW Transparent 6 Gallon
PET Carboy and Accessories**



**Dynalon[®]
Labware**



NEW Transparent 6 Gallon PET Carboy and Accessories

When transparency or reduced gas permeability are important, ask for the Dynalon 6 Gallon PET carboy.

Made from 100% pure food grade PET resin, the product contains no plasticizers or BPA. A wide 4 inch mouth provides easy access for cleaning, sanitizing and filling. Molded graduations mark the 5 and 6 gallon levels.

This carboy is available with and without a PP spigot making it convenient for storage and dispensing of aqueous solutions or for fermentation. As there are no ribs in the design, the carboy has smooth internal and external sides that prevent yeast or bacteria from adhering to the product walls. Carboy is taste and odor free, stain resistant and has excellent gas barrier properties that help prevent oxidation and spoilage.

Each case of 1 each includes instructions on cleaning and sanitizing along with several accessories for use in fermentation applications. These are:

- two polypropylene screw closures, one with a hole that accepts a #10 stopper
- one No. 10 silicone rubber stopper with a hole (useful for a 3 piece fermentation lock)
- one EPDM O-ring that can be inserted into the closure for additional security against leakage.

405564	Description	Dimensions (diam x H)
-006L	Carboy Bottle PET 6 Gallon	11 x 19.3"
-006S	Carboy Bottle PET w Spigot 6 Gallon	12.4 x 19.3"
-0010	Screw Closure Solid PP	5.1 x 1.3"
-0011	Screw Closure PP w/ Hole for #10 Stopper	5.1 x 1.3"
-0012	Stopper #10 with Hole Silicone	2 x 1"
-0013	Stopper #10 Solid Silicone	2 x 1"
-0014	Spigot PP 3/8 inch	3.6 x 1.4 x 2.75"
-0015	O-Ring EPDM, use with Screw Closure	4.4 x 0.1"



405564-0010



405564-0012



405564-0014



405564-0011



405564-0013




405564-0015



405564-006L



405564-006S



**Dynalon
Labware**

175 Humboldt Street, Suite 300
Rochester, NY 14610
Telephone 800-334-7585
Fax 585-334-0241
www.dynalon.com