


Bulk and Gross Packs Snap-Cap Containers



Dynalon®
Labware



Rigid-wall styrene containers with tight fitting polyethylene snap-on caps. Excellent for storing dry materials, powders, tablets, etc. Approximately 16 drams per ounce. 7 dram size may be used with automatic cell counters for collection of sample and, where venous samples are used, for dilution purposes. Max temp PS is 70°C.

426364	Capacity drams (mL)	Height w/o Cap x ID mm	Case Qty
-03	3 (11)	48 x 17	144
-05	5 (18)	50 x 22	144
-07	7 (26)	50 x 26	144
-10	10 (37)	82 x 27	144
-13	13 (48)	82 x 27	144
-15	15 (57)	62 x 33	144
-20	20 (74)	63 x 38	144
-30	30 (111)	75 x 46	72
-40	40 (148)	98 x 45	72

Bulk Vials	Capacity drams (mL)	Bulk Qty	Bulk Price \$	Bulk Caps	Bulk Qty
7003001	3 (11)	2,250	360.00	8003001	14,000
7005011	5 (18)	1,400	233.34	8003010	10,000
7007010	7 (26)	1,000	193.33	8003020	7,500
7010020	10 (37)	600	160.00	8003020	7,500
7013020	13 (48)	500	125.00	8003040	6,000
7015001	15 (57)	500	131.67	8003050	5,000
7016020	20 (74)	400	147.00	8003060	3,000
7030001	30 (111)	230	98.90	8003070	2,500
7040001	40 (148)	175	105.05	8003070	2,500